ES

F313A  
Manual de instrucciones  
WE-2018-F313A

# GALVÁNICO

## ADVERTENCIA

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# GALVANIC

## WARNING

• Never under any circumstances attempt to open or inspect the internal components or accessories of the machine. If an inspection becomes necessary, please contact your supplier or a suitably competent technician.

• Never use the unit in humid or poorly ventilated environments, such as wet rooms, bathrooms etc. to prevent damage to the electrical components of your unit.

• Never use the unit outdoors; exposure to the elements can cause damage.

• If the power cord was destroyed, please ask a professional to service and replace it to avoid hazards, or contact the dealer for service.

## CONTRAINDICATIONS

• Pacemaker (any electrical or rhythmic heart-regulating device)

• Pregnancy

• Thrombosis/Phlebitis

• Skin diseases

• Operations (within the last 6 months)

• Silicone implants

• Hypersensitive or vascular skin

## EC-DECLARATION OF COMPLIANCE

This product meets the following requirements:

1. Conditions set out in the Low Voltage Directive: 2014/35/EU

2. Conditions set out in the Electromagnetic Compatibility Directive: 2014/30/EU

## RECYCLING

 Whenever any of our electrical products or their presentation packaging bear this symbol, it means that these products should not be treated as conventional household waste in Europe. To ensure proper treatment of this waste,

dispose of it in accordance with local laws or as required for the disposal of electrical equipment. This initiative will help preserve save natural resources and improve standards of environmental protection in regard to the treatment of electrical waste.

## WELCOME!

Thank you very much for purchase. Please read the instructions carefully before operating the instrument. You can also ask a professional for guidance to make sure you are using the instrument correctly. We will accept no responsibility for any accidents that occur due to incorrect operation.






The company reserves the right to amend the details of the instruments with no prior notice. If any mistakes are found in these instructions, please feel free to send us an e-mail informing us so that we can correct them.

## PACKAGE CONTENTS


• These components are subject to change without prior notice.

• Appearance of product may differ from the image.

No.	PRODUCT	REF.	QTY.
1	GALVANIC HANDLE (ROLLER)	10042	2
2	GALVANIC HANDLE (BALL)	10043	2
3	GALVANIC WIRES	-	1
4	MAIN UNIT	-	1
5	FEEDER CABLE	-	1

## FRONT PANEL

## INSTRUCTIONS OF USE

1. Connect the power supply via the rear panel.

2. Connect the accessories properly in the lateral panel.

### BEFORE TREATMENT:

1. The treatment area must be washed and dried before any products can penetrate into the tissue.

2. Please make sure that you are working with both hand pieces in contact with the skin of the patient so that the circuit gets closed and the galvanic treatment is effective.

3. The polarity of the probes can be changed by adjusting the "LEAD" button between (+R) or (-G).

4. Cover the roller/ball with a piece of cotton or a sponge that has been dipped into a suitable lotion (according to the treatment). The lotion can also be applied directly to the skin.

### STARTING TREATMENT:

1. Choose the desired galvanic probes (roller or ball). Then properly connect two identical probes in the wires.

2. Turn the unit on by pressing the "POWER" button.

3. Select the desired power output level by pressing "ENERGY" and the desired time by pressing "TIME". You can change any of these settings during the treatment.

4. Press the "START/PAUSE" button.

## CAUTIONS

• Always make sure customers remove any jewelry, and confirm that they have not had any metal implants (pins, joints, etc.) installed since their last treatment.

• Pregnant women and people with sensitive, injured, inflamed skin or blood capillary expansions should not receive this treatment.

• Customers with health issues such as heart disease or infectious diseases and hyper-allergic or hyper-sensitive skin (see contra-indications) may not use this instrument.

• Always make sure the skin surface has been moistened with plenty of conductive product.

## MANINTENANCE

• Take care not to damage the probe rollers. Never use abrasive substances to clean the probe rollers, since damage to the surface of the roller may scratch the surface of the skin during treatment; any damaged areas may also harbor bacteria and make sterilization more difficult.

• After treatment, all metal roller parts must be cleaned, sterilized and dried. Then store them in the accessory holder to avoid any damage.

## ACCESSORIES

REF.	DESCRIPTION
10042/10043	GALVANIC HANDLE
	GALVANIC WIRES

## FAQ

Check the following if the rollers on the probe fail to rotate.

• Dismantle and clean thoroughly.

• If the roller still fails to rotate freely, contact your sales agent.

## TECHNICAL SHEET

Weights and sizes	Packing size (cm)	37x27,5x14
	Gross weight (kg)	2,7
	Product size (cm)	20,5x22x7
Working characteristics	Preset working time	15 min
	Working power	5 W
	Working intensity	30 - 650 micro Amperes
Electrical characteristics	Rated voltage	220V ~ 240V
	Frecuency	50Hz/60Hz
	Rated power	29 W

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