

ORGANIC

Elegant stool ergonomically designed with a saddle-shaped seat to ensure maximum comfort. Equipped with a gas lift to adjust height. White PU upholstery of high quality and very stable chrome base.

Taburete de diseño ergonómico y elegante con asiento en forma de silla de montar. Altura regulable mediante pistón de gas. Tapizado blanco en PU de alta calidad y base cromada.

Tabouret de design ergonomique et élégant. Siège en forme de selle de cheval. Hauteur réglable par vérin à gaz. Revêtement en PU blanc de haute qualité et base chromée.

Sattelhocker mit einem ergonomischen und eleganten Design. Die Höhe kann mittels Gasdruckfeder reguliert werden. Hochwertige PU-Polsterung und Chromfuß.

Sgabello dal design ergonomico ed elegante e sedile a forma di sella. Regolabile in altezza con un pistone a gas. Rivestimento in PU bianco di alta qualità e stativo cromato.

Saddle-shaped Stool
Ref. 1022A
Weight: 6 Kg
Upholstery: PU

1 Up and Down 59/78 cm
2 Base diameter 46 cm

OPTIONS
Lower gas lift Ref. 30004
Up and Down 55/73 cm

COLOURS



PACKING SIZE cm
48x22x48
N.W. kg
6
G.W. kg
7

