

## CLAVI

**3-section chair with 3 motors that control height, backrest and legrest inclination. Reclining armrests and adjustable headrest with breather hole. Also available the 4-motor model (WKE015.4.A26) that allow a motorised adjustment of the seat inclination. White PVC upholstery that eases cleaning. Roll holder included.**

Sillón de estética de 3 cuerpos con 3 motores que controlan la altura, la inclinación del respaldo y del reposapiés. Reposabrazos abatibles y apoyacabezas regulable con agujero facial. Disponible también con 4 motores (modelo WKE015.4.A26) que permiten la regulación motorizada del asiento. Tapizado blanco en PVC y fácil de limpiar. Portarrollos incluido.

Fauteuil à 3 plans équipé de 3 moteurs qui contrôlent la hauteur et l'inclinaison du dossier et du repose-jambes. Accoudoirs rabattables et repose-tête réglable avec trou visage. Il est disponible sur 4 moteurs (modèle WKE015.4.A26) pour le réglage motorisé de l'inclinaison du siège. Revêtement en PVC blanc et facile à nettoyer. Porte-rouleau inclus.

*Dreiteilige Liege. Mit 3 Motoren um die Höhe, die Neigung der Rückenlehne und des Fußteiles zu regulieren. Armlehnen beidseitig abklappbar und regulierbare Kopfstütze mit Gesichtsausschnitt. Auch mit 4 Motoren erhältlich (Modell WKE015.4.A26), um die Neigung des Sitzes zu regulieren. Weiße PVC Polsterung, leicht abwaschbar. Mit Papierrollenhalter.*

*Poltrona estetica a 3 corpi con 3 motori che controllano l'altezza, l'inclinazione dello schienale e del poggiatesta. Braccioli pieghevoli e poggiatesta regolabile con foro facciale. Disponibile anche con 4 motori (modello WKE015.4.A26) che consentono la regolazione motorizzata del sedile. Rivestimento in PVC bianco facile da pulire. Portarotolo incluso.*

### Electric Beauty Bed

Ref. WKE014.3.A26

Size without armrest: 203x63x65/93 cm

Size with armrest: 203x84x65/93 cm

Weight: 72 Kg

Upholstery: PVC

Motors: 3

### Electric Beauty Bed

Ref. WKE015.4.A26

Size without armrest: 203x63x65/100 cm

Size with armrest: 203x84x65/100 cm

Weight: 74 Kg

Upholstery: PVC

Motors: 4

**1** Breather hole with headrest cushion

**2** Angle 0°-68°

**3** Angle -10° (Only 4 motors)

**4** Angle 0°-50°

**5** Extension ±17 cm

**6** Up and Down 65/93 cm

**7** Paper roll holder

**8** Hand controller

### OPTIONS

**Bed cover**

**9** Foot controller Ref. FC003.WK/FC004.WK

### PACKING SIZE cm

163x68x70,5

N.W. kg

72/74

G.W. kg

87/89

MAX WEIGHT

180 kg

NEW

3 - 4  
MOTORS

7



8



4



1

5

2

3

6

9



FC003.WK



# CLAVI

